

## **Drill 1-Dexterity Skills**

### **Drill Requirements:**

Abdomen Model 1

### **Instruments:**

Atraumatic bowel grasper

Maryland Dissector

### **Other Equipment:**

Video monitor

Camera

Telescope 4mm, 30 degree

### **Task A: Ring to Vertical Peg Transfer**

**Objective:** Transfer 8 rings onto a vertical peg

### **Performance Scoring:**

Time to complete task

Total number of errors

    Dropped ring, missed attempt

    Instrument withdrawn/removed from trainer

### **Task B: Bead to String Transfer**

**Objective:** Transfer 8 beads onto a string, alternating colors

### **Performance Scoring:**

Time to complete task

Total number of errors

    Dropped bead, missed attempt

    Instrument withdrawn/removed from trainer

    Bead colors do not alternate

### **General Principles:**

Time starts when participant places instruments into trainer, instruments are visible, and participant expresses that they are ready to begin. Each task is timed individually and start of timing is same for each task.

Errors are recorded if the participant drops a ring or bead or if they release ring or bead prior to it being placed on peg or string or if instruments are inadvertently pulled from the box trainer (tip of instrument out of box—one error for each instrument removed i.e., both instruments pulled out = 2 errors). In the bead task, failing to alternate bead colors is scored as an error each time the bead color is incorrect.

Time to complete each task as well as the number of errors that occurred during completion of each task are recorded.

# Neonatal MIS Trainer Project Participant Score Sheet

Session Date: \_\_\_\_\_ REDCAP Record ID: \_\_\_\_\_ - \_\_\_\_\_

Session Type (select one):

Evaluation Session (TS-E)				
<input type="checkbox"/> Baseline	<input type="checkbox"/> 6-month	<input type="checkbox"/> 12-month	<input type="checkbox"/> 18-month	<input type="checkbox"/> 24-month
<input type="checkbox"/> Mandatory Practice Session (TS-P)		<input type="checkbox"/> Optional Practice Session (TS-O)		

## Drill 1-Dexterity Skills

### Task A: Ring to Peg Transfer

Time to transfer 8 rings onto a vertical peg

\_\_\_\_\_ Minutes      \_\_\_\_\_ Seconds

**Errors\***

Dropped ring	Count:	
Missed attempt	Count:	
Instrument(s) pulled out of trainer	Count:	

### Task B: Bead to String Transfer

Time to transfer 8 beads onto a string, alternating colors

\_\_\_\_\_ Minutes      \_\_\_\_\_ Seconds

**Errors\***

Dropped bead	Count:	
Missed attempt	Count:	
Instrument(s) pulled out of trainer	Count:	
Bead color does not alternate	Count:	

\*Record each error and indicate **total** number of each error type